



**CERTIFICATE OF ANALYSIS**  
**Phosphate - Buffered Saline Solution**  
**(Catalog Number 600-108)**

Lot Number: M122180

Date of Manufacture: Dec2022

Product Expiry: Dec2024

Storage Temperature: Ambient

**Product description:** *Phosphate - Buffered Saline Solution is sterile-filtered through a 0.2 µm filter prior to storage.*

Test	Methodology	Specification	Analysis
<b>Biological Testing</b>			
Endotoxin	USP <85>	Report	<0.750 EU/mL
<b>Microbiological Testing</b>			
Sterility	USP <71>, EP 2.6.1		
Bacteria		No Growth	No Growth
Fungi		No Growth	No Growth
<b>Physical Testing</b>			
pH@RT	USP<791>	Report	7.45
Osmolality	USP<785>	Report	281 mOsm/kg

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

VANESSA YOUNG  
Name

QA ASSOCIATE II  
Title

06 JAN 2023  
Date