



CERTIFICATE OF ANALYSIS
Polymixin B Solution (50mg/mL)
(Catalog Number 200-283-0100mL)

Lot Number: F021070

Date of Manufacture: Oct2021

Product Expiry: Oct2023

Origin: United States

Storage Temperature: 2-8°C

For Cell culture, Further Manufacturing or Research use Only. Not for Direct Therapeutic Use.

Product description: Polymixin B Solution is sterile-filtered through a 0.2 µm filter prior to refrigerating.

Test	Methodology	Specification	Analysis
Microbiological Testing			
Sterility	USP<71>, EP 2.6.1		
Bacteria		No Growth	No Growth
Fungi		No Growth	No Growth
Physical Testing			
Osmolality	USP<785>, EP 2.2.35	Test and Report	53 mOsm/kg
pH	USP<791>	Test and Report	5.89
Active Unit Activity			
Polymixin B Sulfate Powder			7925 U/mg

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Hina Z Khan
Name

QC Associate I
Title

13 NOV 2021
Date