



CERTIFICATE OF ANALYSIS

Stasis™ Stem Cell Qualified U.S. Fetal Bovine Serum **(Catalog Number 100-525)**

Lot Number: A424061

Date of Manufacture: Oct2023

Product Expiry: Oct2028

Origin: United States

Storage Temperature: ≤ -10°C

For cell culture, further manufacturing or research use only. Not for direct therapeutic use.

Product Description: Stasis™ Stem Cell Qualified U.S. Fetal Bovine Serum is sterile-filtered through a 0.1 µm filter prior to freezing.

Test	Methodology	Specification	Analysis
Biological Testing			
Endotoxin	USP<85>/Ph. Eur. 2.6.14	<10.0 EU/mL	<0.50 EU/mL
Hemoglobin	USP<857>	<20.0 mg/dL	4 mg/dL
Microbiological Testing			
Sterility	USP<71>/Ph. Eur. 2.6.1		
Bacteria		No Growth	No Growth
Fungi		No Growth	No Growth
Mycoplasma	USP<63>/Ph. Eur. 2.6.7	Not Detected	Not Detected
Viral Testing	9CFR113.53c [113.46,113.47]		
Cytopathic Agents		Not Detected	Not Detected
Hemadsorbing Agents		Not Detected	Not Detected
Extraneous Viruses (FAB)			
Bluetongue Virus		Not Detected	Not Detected
Bovine Adenovirus		Not Detected	Not Detected
Bovine Parvovirus		Not Detected	Not Detected
Bovine Viral Diarrhea (Cytopathic)		Not Detected	Not Detected
Bovine Respiratory Syncytial Virus		Not Detected	Not Detected
Rabies Virus		Not Detected	Not Detected
Reovirus		Not Detected	Not Detected
Physical Testing			
Electrophoretic Pattern		Normal	Normal
Osmolality	USP<785>/Ph. Eur. 2.2.35	260 – 350 mOsm/kg	311 mOsm/kg
pH	USP<791>/Ph. Eur. 2.2.3	6.5 – 8.5	7.38
Biochemistry Testing			
Albumin		Test and Report	2.2 g/dL
Alkaline Phosphatase		Test and Report	256 U/L
ALT (SGPT)		Test and Report	3 U/L
AST (SGOT)		Test and Report	20 U/L
GGT		<10 U/L	4.0 U/L
Bilirubin, Total		Test and Report	0.1 mg/dL
BUN		Test and Report	14 mg/dL



CERTIFICATE OF ANALYSIS

Stasis™ Stem Cell Qualified U.S. Fetal Bovine Serum

(Catalog Number 100-525)

Lot Number: A424061

Test	Methodology	Specification	Analysis
Calcium		Test and Report	12.7 mg/dL
Chloride		Test and Report	97 mmol/L
Cholesterol		Test and Report	32 mg/dL
Cortisol		Test and Report	<0.20 µg/dL
Creatinine		Test and Report	2.30 mg/dL
Glucose		Test and Report	91 mg/dL
IgG		<0.300 mg/mL*	0.1 mg/mL
Insulin		Test and Report	4.23 µIU/mL
Iron, Serum		Test and Report	169 µg/dL
Phosphorus		Test and Report	10 mg/dL
Potassium		Test and Report	13.1 mmol/L
Progesterone		Test and Report	<0.06 ng/mL
Protein, Total	USP<1057> Ph. Eur. 2.5.33	2.5 – 6.0 g/dL	3.2 g/dL
Sodium		Test and Report	134 mmol/L
Testosterone		Test and Report	0.02 ng/mL
T3		Test and Report	114 ng/dL
T4		Test and Report	10.7 µg/dL
Triglycerides		Test and Report	80 mg/dL
Uric Acid		Test and Report	2.1 mg/dL

Our fetal bovine serum meets all USDA requirements and has passed ante and post-mortem inspection. It is collected in countries recognized as being free of foot-and-mouth disease and rinderpest. The origin of all fetal bovine serum is traceable by lot number, date, and country of collection. Precipitates may develop in this product upon freezing and/or thawing; this occurrence does not impact culture performance. *Some geographic regions may test higher than 0.300 mg/mL. If this be the case, refer to GGT results less than 10 U/L to ensure FBS purity.

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and has confirmed that the testing meets the specifications presented on this Certificate of Analysis.

Karen Zuniga
Name

25NOV2024
Date

QA Associate II
Title