



CERTIFICATE OF ANALYSIS

Yssel's Custom T-Cell Medium without Phenol Red (Catalog Number 900-133)

Lot Number: M425004

Date of Manufacture: Feb2025

Storage Temperature: $\leq -10^{\circ}\text{C}$

Product Expiry: Feb2027

For cell culture, further manufacturing or research use only. Not for direct therapeutic use.

Test	Methodology	Specification	Analysis
Biological Testing			
Endotoxin	USP<85>, EP 2.6.14	<5 EU/mL	<1.00 EU/mL
Microbiological Testing			
Sterility	USP<71>, EP 2.6.1		
Bacteria		No Growth	No Growth
Fungi		No Growth	No Growth
Physical Testing			
pH	USP<791>	6.5 – 7.5	7.01
Osmolality	USP<785>, EP 2.2.35	250 – 350 mOsm/kg	294 mOsm/kg

The testing that has been performed as part of this lot release has been reviewed by Quality Assurance personnel and it has been confirmed that the testing meets the specifications presented on this Certificate of Analysis. This product was manufactured at 3230 Reed Ave, West Sacramento, CA, U.S.A., which is an animal origin facility. All final product processes/packaging were also performed at 3230 Reed Ave, West Sacramento, CA, U.S.A., an animal origin facility. This facility is ISO 13485:2016 certified.

Korrina Nicewonger
Name

26 FEB 2025
Date

QA Associate II
Title